

# 90101212 SIMPLE 3

## FUNCTIONS



Backrest adjustment : 0 - 70°



Thigh rest adjustment : 0 - 40°



Foot rest mechanical adjustment : 0 - 25°



Vascular position



Fowler position



Trendelenburg position : 0 - 12°


# HOSPITAL BED WITH 2 CRANKS

## OPTIONS


### GROUP 1 / CASTOR OPTIONS



### GROUP 2 / SIDE RAIL OPTIONS



### GROUP 3 / OTHER ACCESSORY OPTIONS



## DESCRIPTION OF OPTION CODES

<b>T1</b>	TWIN CASTOR TWO WITH BRAKE (120 MM)	<b>W0</b>	WITHOUT OTHER ACCESSORIES
<b>T16</b>	CASTOR TWO WITH BRAKE (125 MM)	<b>W2</b>	LIFTING POLE : BEARING CAPACITY 75 Kg
<b>T4</b>	METAL CASTOR TWO WITH BRAKE (125 MM)	<b>W4</b>	TELEVISION STAND FOR 17" LCD TV
<b>K0</b>	WITHOUT SIDE RAIL	<b>W5</b>	TRACTION APPARATUS
<b>K1</b>	WITH METAL SIDE RAIL	<b>W8</b>	IV STAND : ONE HEIGHT ADJUSTMENT SYSTEM
<b>K3</b>	WITH METAL SIDE RAIL		

Must select only one of them.

Do not select the others if you select without accessory. If you select the others you can select all or 1-2-3 ... of them.

## GENERAL SPECIFICATIONS

- Easily cleanable metal bed platform designed for the ventilation of the bed.
- Head and foot movements with rotating levers.
- Trendelenburg movement by the help of shock absorber mechanism at foot side.
- Ergonomic, durable and easily cleanable metal rails system with maximum safety.
- Melamine coated chipboard head and foot panels.
- Removable and installable head system that can be also used as a resulation board if necessary.
- Shock resistance plastic bumper wheels.
- Metal components dyed with anti-corrosion electrostatic powder paint after being treated with phosphate.
- Optional melamine coated chipboard color choice.<sup>(RSX)</sup>

## TECHNICAL SPECIFICATIONS

■ External dimension	:	105 x 204 x 50 cm
■ Recommended mattress dimension	:	85 x 190 x 12 cm
■ Safe Operating Load (Inc. Mattress and Accessories)	:	230 kg
■ Weight	:	~ 82 kg
■ Volume	:	1,17 m <sup>3</sup>



\* RS : Melamine color options page 81.